

# Molecular approaches: TB day

Monday 2 September

# MOLECULAR ASSAYS FOR TUBERCULOSIS: increased levels of complexity

INFORMATION

*Mtb*/no *Mtb*  
RIF

*Mtb*/no *Mtb*  
15 drugs  
Genotype

*Mtb*/no *Mtb*  
all drugs  
Genotype  
Transmission

TARGETS

MTB-specific  
1 (4) drug-target(s)

MTB-specific  
18 drug-targets

4.4Mb  
All  
drug-targets

SAMPLE  
TYPE





Raw  
specimen

Treated  
Specimen  
\*\*\*

Cultured  
isolate

MOLECULAR  
ASSAY

Simple  
NAATs




# Cepheid GeneXpert® System

Nucleic acid amplification tests (NAATs)



# The impact of Xpert on TB diagnosis and care

Impact on Patient Isolation Pathway  
"Test and Release" or "Test and Treat"



# System vs. Assays - semantics

- **THE SYSTEM: GeneXpert® System**
  - “The system is designed to purify, concentrate, detect and identify targeted nucleic acid sequences thereby delivering answers directly from unprocessed samples.”
- **THE ASSAY: Xpert <<your favourite assay>>**
  - Xpert® MTB/RIF
  - Xpert® MTB/XDR
  - Xpert® Xpress Flu
  - Xpert® Xpress Strep A

# GeneXpert systems for different throughput needs

Available in two, four, or 16-module configurations, each system has our proven GeneXpert module at its analytic core and uses patented cartridge technology for every Xpert® test.

Look for the blue band signifying 10-color technology capabilities.



GeneXpert II

GeneXpert IV

GeneXpert XVI

80



Up to 2000 tests/day

GeneXpert Infinity

# How does it work?



# How does the PCR detection used by Xpert MTB/RIF Ultra work?

- **TB or no TB? (That is the question)**
  - Targets:
- **Xpert MTB/RIF:**
  - RIF susceptible or RIF resistant?
    - Target: *rpoB* – RNA polymerase
    - Rifampicin resistance determining region – 82 bp genomic region in *rpoB*
    - Wild type and mutant probes
- **Xpert MTB/XDR**
  - 10 colour-assay
  - Extended drug resistance:
    - INH
    - FQs
    - SLIDs
    - ETH



# Xpert assays for TB



2006

Cepheid, FIND, Rutgers University, and NIAID begin development of **Xpert® MTB/RIF** with funding from the Bill & Melinda Gates Foundation



2009

**Xpert MTB/RIF** launches

2017

**Xpert MTB/RIF Ultra** launches and gets WHO endorsement



2020

**Xpert MTB/XDR** launched on new 10-color multiplexing technology

2021

**Xpert MTB/XDR** gets WHO endorsement

# Interpreting results

110008945 ITG Antwerpen

12/23/22 11:17:22

## Test Report

Patient ID:  
Sample ID: CT202201273  
Test Type: Specimen  
Sample Type:

### Assay Information

Assay	Assay Version	Assay Type
Xpert MTB-XDR	1	In Vitro Diagnostic

Test Result:

MTB DETECTED;

INH Resistance NOT DETECTED;

FLQ Resistance NOT DETECTED;

AMK Resistance NOT DETECTED;

KAN Resistance NOT DETECTED;

CAP Resistance NOT DETECTED;

ETH Resistance NOT DETECTED

# Interpreting results

110008945 ITG Antwerpen

12/23/22 11:18:13

## Test Report

Patient ID:  
Sample ID: CT202201279  
Test Type: Specimen  
Sample Type:

### Assay Information


Assay	Assay Version	Assay Type
Xpert MTB-XDR	1	In Vitro Diagnostic

Test Result:


MTB DETECTED;  
INH Resistance NOT DETECTED;  
FLQ Resistance DETECTED;  
AMK Resistance NOT DETECTED;  
KAN Resistance NOT DETECTED;  
CAP Resistance NOT DETECTED;  
ETH Resistance NOT DETECTED

# PRACTICAL

- Work in groups of 2
- Eyeball sample volume
- Add 2x volume sample reagent
- Mix well
- Incubate at RT for 15 min
- Mix again at 10 min mark
- Add 2 ml to cartridge
- Load into machine
- Follow software prompts
- Wait for assay to run
- Interpret results




# New assays: testing the limits





# New assays: testing the limits

- **Xpert MTB/XDR:**  
Determining the limit of  
FQ heteroresistance detection



# Deeplex Myc-TB

Targeted Next-Generation Sequencing




# tNGS & WGS

Two NGS approaches for DR-TB control


**tNGS**

Focus on genomic  
regions of interest



**WGS**

Comprehensive analysis  
of the entire genome





# tNGS & WGS

Two NGS approaches for characterisation of DR-TB

**tNGS**

Focus on genomic  
regions of interest



**WGS**

Comprehensive analysis  
of the entire genome



# tNGS & WGS: Head to head

Characteristic	tNGS - <i>Deeplex® Myc-TB</i>	WGS
Input sample	Sputum	Culture
Batching	Min 20	Min 4 (Illumina), Min 1 (ONT)

# tNGS & WGS: Head to head


Characteristic	tNGS - <i>Deeplex® Myc-TB</i>	WGS
Input sample	Sputum	Culture
Batching	Min 20	Min 4 (Illumina), Min 1 (ONT)
bp investigated	10,000 - 13000 bp	~ 4 million bp
Drug resistance	12 – 15 drugs fixed	23 drugs flexible
Limit of variant detection	1-3%	1-3% if 1000x depth

# tNGS & WGS: Head to head

Characteristic	tNGS - <i>Deeplex® Myc-TB</i>	WGS
<b>Input sample</b>	Sputum	Culture
<b>Batching</b>	Min 20	Min 4
<b>bp investigated</b>	10,000 - 13000 bp	~ 4 million bp
<b>Drug resistance</b>	12 – 15 drugs fixed	23 drugs flexible
<b>Limit of variant detection</b>	1-3%	1-3% if 1000x depth
<b>Phylogenetic Typing</b>	Moderate resolution (SNPs/spoligotyping)	High resolution
<b>Transmission Inference</b>	Limited	Accurate


# Implications of a culture step

**tNGS - *Deeplex*® *Myc-TB***



# Implications of a culture step

WGS



- Diversity loss
- Increased turnaround times
- Biosafety and infrastructural needs

# A comprehensive approach: Drug resistance

RIF	INH	EMB	PZA	BDQ	LZD	MXF	LFX	STM	AMK	KAN	CAP	CFZ	ETO
<i>rpoB</i>	<i>inhA</i>	<i>embB</i>	<i>pncA</i>	<i>mmpR5</i>	<i>rpIC</i>	<i>gyrA</i>	<i>gyrA</i>	<i>rrs</i>	<i>rrs</i>	<i>rrs</i>	<i>rrs</i>	<i>mmpR5</i>	<i>inhA</i>
	<i>ahpA</i>				<i>rrl</i>	<i>gyrB</i>	<i>gyrB</i>	<i>rspl</i>		<i>eis</i>	<i>tlyA</i>		<i>ethA</i>
	<i>katG</i>							<i>gid</i>					

- 13-15 drugs
- 18 targets



# A comprehensive approach: Drug resistance

RIF	INH	EMB	PZA	BDQ	LZD	MXF	LFX	STM	AMK	KAN	CAP	CFZ	ETO
<i>rpoB</i>	<i>inhA</i>	<i>embB</i>	<i>pncA</i>	<i>mmpR5</i>	<i>rplC</i>	<i>gyrA</i>	<i>gyrA</i>	<i>rrs</i>	<i>rrs</i>	<i>rrs</i>	<i>rrs</i>	<i>mmpR5</i>	<i>inhA</i>
<i>rpoA</i>	<i>ahpA</i>	<i>embA</i>	<i>clpC1</i>	<i>mmpS5</i>	<i>rrl</i>	<i>gyrB</i>	<i>gyrB</i>	<i>rspl</i>	<i>whiB7</i>	<i>eis</i>	<i>tlyA</i>	<i>mmpS5</i>	<i>ethA</i>
<i>rpoC</i>	<i>katG</i>	<i>embC</i>	<i>panD</i>	<i>mmpL5</i>				<i>gid</i>	<i>eis</i>	<i>whiB7</i>	<i>whiB6</i>	<i>mmpL5</i>	<i>ethR</i>
<i>Rv2752c</i>	<i>mshA</i>	<i>embR</i>	<i>Rv1258c</i>	<i>atpE</i>				<i>whiB7</i>	<i>whiB6</i>		<i>ccsA</i>	<i>atpE</i>	<i>mshA</i>
	<i>Rv1258c</i>	<i>ubiA</i>	<i>PPE35</i>	<i>pepQ</i>				<i>Rv1258c</i>	<i>ccsA</i>		<i>fprA</i>	<i>pepQ</i>	<i>Rv3083</i>
	<i>ndh</i>		<i>Rv3236c</i>	<i>Rv1979c</i>				<i>whiB6</i>	<i>fprA</i>		<i>aftB</i>	<i>Rv1979c</i>	<i>ndh</i>
	<i>Rv2752c</i>								<i>aftB</i>				

## DLM

*fdg1*

*ddn*

*fbiA*

*fbiB*

*fbiC*


*fbiD*



World Health Organization




# PRACTICAL – Deeplex Myc-TB



# Deeplex® Myc-TB workflow

## PRACTICAL

## LAB MANUAL: P19



# Oxford Nanopore Technologies

Whole genome sequencing






Please make an account on Galaxy before Thursday  
afternoon's session

<https://usegalaxy.org/>

# How ONT sequencing works - recap



# Overview of ONT sequencing approaches

## Investigations

---

Structural variation

Gene expression

Splice variation

Fusion transcripts

Single cell

SNVs and phasing

Identification

Assembly

Epigenetics

Chromatin  
conformation

## Techniques

---

Whole genome

Targeted

Whole transcriptome

Metagenomics

Short fragment mode

# ONT solutions for WGS

Application-free, native DNA sequencing and retained base modification			Amplification-based for low DNA amounts or quality	
	Ligation Sequencing Kit	Rapid Sequencing Kit	Ultra-Long Sequencing Kit	Rapid PCR Barcoding Kit
<b>Preparation time</b>	60 min	10 min	200 min + O/N elution	15 min + PCR
<b>Input requirement</b>	1,000 ng gDNA; 100–200 fmol amplicons or cDNA	200 ng gDNA (<30 kb)	6 M cells	1–5 ng gDNA
<b>Fragmentation</b>	Optional	Transposase based	Transposase based	Transposase based
<b>Read length</b>	Equal to fragment length	Random distribution, dependent on input fragment length	50–100+ kb N50	~2 kb
<b>Typical throughput</b>	● ● ●	● ● ○	● ● ○	● ● ○
<b>Multiplexing options</b>	24 plex, 96 plex	24 plex, 96 plex	-	12 plex
<b>Methylation included</b>	Yes	Yes	Yes	-
<b>Overview</b>	Optimised for output; retained base modifications; control over read length	Simple and rapid; retained base modifications	Optimised for production of ultra-long reads (N50 ≥50 kb); retained base modifications	Simple and rapid

# Rapid barcoding kit overview


## Tagmentation

Fragmentation and adapter ligation in a single step

## Sample pooling

## Sequencing adapter ligation

Starting material





# Rapid barcoding kit overview

Library preparation step	Process	Time	Stop option
DNA barcoding	Tagmentation of the DNA using the Rapid Barcoding Kit V14	15 minutes	4°C overnight
Sample pooling and clean-up	Pooling of barcoded libraries and AMPure XP Bead clean-up	25 minutes	4°C overnight
Adapter ligation	Attach the sequencing adapters to the DNA ends	5 minutes	We strongly recommend sequencing your library as soon as it is adapted
Priming and loading the flow cell	Prime the flow cell and load the prepared library for sequencing	5 minutes	

# PRACTICAL – stations and sequencing

- **24 samples to be prepared** (200 ng gDNA aliquoted)
- **3 library pools: 3 x stations for library preparation**
  - Each group to pool their samples – 8 samples per pool
  - Each sample = one *Mycobacterium tuberculosis* clinical isolate
- **3 x libraries will be loaded for sequencing – training room**
  - 3 x Laptops set up with MinKNOW software and MinION sequencers
  - 3 volunteers to load flowcells

# PRACTICAL – flowcell check and loading (seminar room)

- **Quality control of flow cells – flowcell health (done)**
  - Check how many active pores:
    - Up to 2048 pores
    - Ideally more than 800 active pores
- **Follow in instructions in lab manual: p25 - 26**
  - Prime flowcell
  - Prepare library for loading for loading
  - Load the flow cell and start the sequencing


# PRACTICAL – priming and loading the flow cell



# Tips and tricks for loading a flow cell

- **Removing liquid from the priming port to ensure that there are no bubbles in the sequencing channel**
  - Insert the tip of a P1000
  - Slowly turn the dial to remove  $\pm 200\mu\text{l}$  liquid
- **Loading the flowcell: work quickly and accurately**
  - Loading beads precipitate quickly
  - The order in which to open/close the ports matter
- **Reverse pipetting**
  - Use RP to add all reagents to the flow cell

# Reverse pipetting



**Collect reagent from reagent tube**

**Load only indicated volume into flow cell**

**Empty remaining volume left in tip  
in reagent tube**