



University of Antwerp
| GHI | Global Health Institute

Next Generation Sequencing Technologies

Anzaan Dippenaar

31 August 2024




DNA sequencing

A BRIEF HISTORY



Friedrich Miescher:

- isolated various phosphate-rich chemicals from white blood cells, which he called nuclein



Avery–MacLeod–McCarty:

- DNA is the “genetic material”




Watson & Crick:

- Double Helix Structure of DNA


History of DNA sequencing

1977




Maxam-Gilbert:
DNA sequencing by
chemical
degradation

Frederick Sanger *et al*:
DNA sequencing with
chain-terminating
inhibitors




1981



Messing *et al*:
shotgun sequencing of
cauliflower mosaic virus

1983



Kary Mullis:
Polymerase Chain
reaction (PCR)


The history of sequencing – 1st gen

1987




First automated sequencer: ABI Prism 370

1995



Craig Venter and Hamilton Smith:
first complete genome of a
free-living organism: *H.*
influenza

2001



Draft of the human
genome published

The history of sequencing: NGS

2005




Roche 454
pyrosequencing

2006



Solexa (Illumina)
Next Generation
Sequencer

2007



Applied Biosystems
Next Generation
Sequencer : SOLiD

2011



Ion Torrent: Personal
Genome Machine
(PGM)

The history of sequencing: 3rd gen

2009



Helicos single molecule
sequencer : Helicos
Genetic Analyser System

2011



Pacific Biosciences single
molecule sequencer : PacBio RS
Systems

2012



Oxford Nanopore
Technologies -
The GridION and MinION

DNA sequencing

THE PRINCIPLES

Next-Generation Sequencing platforms

- **Different approaches are used:**
 - 1) Sequencing by synthesis - uses DNA polymerase (e.g. Illumina, Pyrosequencing)
 - 2) Sequencing by ligation - uses DNA ligase (e.g. SOLiD)
 - 3) Direct molecule sequencing (e.g. ONT)

Choosing the most appropriate platform


Platform	Primary advantages	Primary disadvantages
Sanger Sequencing	Low cost for small studies	High cost for large amounts of data
Illumina	Moderate run cost, low error rate, various read length configurations, standardised data format	High instrument costs, relatively short reads, large amounts of data produced*
PacBio	Longest available read length (single molecule real time sequencing), short run time	Expensive equipment and per-sample costs, high(er) error rate
Ion Torrent	Low-cost instruments and disposable chips, simple machine	High error rate, fewer reads at higher cost per Mb relative to e.g. Illumina MiSeq, smaller user community, data format
Oxford Nanopore Technologies	Minlon – small, portable instrument, low capital (hardware) cost, long reads, error rate does not increase along read length	High error rate, limited library preparation approaches, high amount of input DNA, analysis not standardized

Understanding sequencing – PCR



Sequencing by synthesis

1. Fragment genomic DNA
2. Add 2x different synthetic adapters complementary to that on the flowcell.
3. Hybridize to solid flowcell.
4. Perform bridge amplification until flowcell is covered with DNA strands.
5. Target strands will now be clustered on the flowcell.
6. Add adapter specific primers with fluorescently labeled DNA basis (ATGC)/
7. As a base is used during amplification, a signal is detected and captured.
8. These contigs are then compared to a reference genome.



<https://www.youtube.com/watch?v=CZeN-IgjYCo>

<https://www.illumina.com/science/technology/next-generation-sequencing/beginners/ngs-workflow.html#:~:text=The%20next%2Dgeneration%20sequencing%20workflow,to%20plan%20your%20NGS%20workflow.>

Sequencing by synthesis


<https://www.youtube.com/watch?v=CZeN-IgjYCo>

<https://www.illumina.com/science/technology/next-generation-sequencing/beginners/ngs-workflow.html#:~:text=The%20next%2Dgeneration%20sequencing%20workflow,to%20plan%20your%20NGS%20workflow.>

Nanopore DNA Sequencing

Oxford Nanopore's technology uses nanopores embedded in a lipid membrane.


As a DNA molecule moves through the nanopore, minute changes in the electrical current are recorded, and later interpreted.



Nanopore DNA Sequencing


ONLINE VIDEO: ONT

NGS analysis approaches



Reference mapping approach

- Known genome sequence of closely related organism
- Analysis is less resource-intensive
- Limited to known genomic regions



De novo assembly

- Limited knowledge of genome sequence or closely related organism
 - Analysis is resource-intensive
- Describe novel genomic structures
 - Species characterisation

WGS ANALYSIS



Quality control and preprocessing



Alignment to a reference genome



Refine and compress alignment



Variant identification



Annotation




Next generation sequencing
of *M. tuberculosis*

FROM SAMPLE TO SEQUENCE

Microbiology

Clinical specimen




Liquefaction




Decontamination

Mtb Enrichment



Inactivation




Molecular biology

Mtb lysis




DNA extraction and quality control



Library preparation and sequencing




Data analysis




Bioinformatics

NGS APPLICATIONS FOR TB RESEARCH AND CLINICAL CARE

Surveillance




Epidemiology




Drug resistance





Heteroresistance




Transmission



TARGETED vs. WHOLE GENOME SEQUENCING



Isoniazid
Rifampicin
Bedaquiline




Logistics and needs for TDS vs standard WGS vs advanced WGS

	Targeted Sequencing (Deeplex-MycTB)	Standard WGS	Advanced WGS
Application	Culture-free prediction of resistance to anti-TB drugs	Culture-based prediction of resistance to anti-TB drugs	Culture-based prediction of resistance to anti-TB drugs, transmission dynamics
Input material	DNA extracted directly from specimens (if bacterial load high) or cultured isolates	High quality purified DNA extracted from cultured <i>Mtb</i> isolates	High quality purified DNA extracted from cultured <i>Mtb</i> isolates
Staff requirements	Sequencing technician	Sequencing technician	Technician and Bioinformatician
Depth of coverage	> 1000 X (up to 100k X)	>10x	>30x (depending on application)
TAT (assuming in-house sequencer)	2 days	3-40 days culture 1 day DNA extraction 2-3 days sequencing 1 day analysis (Need for culture increases TAT)	3-40 days culture 1 day DNA extraction 2-3 days sequencing 2-day analysis (standard) 4+ days analysis (additional) (Need for culture increases TAT)
Data analysis	User friendly cloud-based system	User-friendly graphical interface	Bioinformatics analysis tools. (Graphical interface, command line)

NGS - Applications in *M. tuberculosis*
tNGS vs WGS

PROJECT OVERVIEW: deepMTB

- Impact of (Hetero)resistance on treatment success of DR-TB




Comparison of
Hetero-resistance detected by different approaches
+

Quantify prevalence culture bias

TARGETED & WHOLE GENOME SEQUENCING

PROJECT OVERVIEW: deepMTB



- Lab-made resistant and susceptible *M. tuberculosis* mixtures
- Available from BE-based bacterial biobank
- Developers of diagnostic tests




NGS - Applications in *M. tuberculosis*

NGS has higher resolution compared
to classic genotyping techniques

WGS: MDR Beijing outbreak in CAPE TOWN


- 1993-2015
- 84 patients (127 isolates)
- Identical Beijing IS6110-RFLP
- Resistant to 4 frontline drugs
- WGS revealed characteristic: *pncA* 100ATC; 160GCA mutations



WGS of clinical primary *Mtb* liquid cultures for

PRECISION MEDICINE

PROJECT OVERVIEW: SMARTT



- **Pragmatic randomised control trial**
 - WGS and Standard of care
- **WGS-based comprehensive drug-resistance profile to guide DR-TB treatment**
 - Individualised treatment / personalised medicine
- **Low-middle income country**
 - Freestate Province, South Africa
- **Primary outcome**
 - Bacteriological response to treatment

PROJECT OVERVIEW: SMARTT



MYCOBACTERIUM TUBERCULOSIS GENOME SEQUENCING REPORT

SMARTT
Sequencing Mycobacteria and Algorithm-based Resistant Tuberculosis Treatment

Patient Name		Phone	
Birth date		Sample	LP250
Gender	MALE	Weight (Kg)	55.0
Facility		Facility contact	
Requested by		Sequenced from	MGIT
Sample type	Sputum	Sample date	XXX
Sample number		Report date	2022-11-01

Drug Resistance Profile

This *Mycobacterium tuberculosis* strain is predicted to be:

Resistant to: **Ethambutol**, **Ethionamide**, **Prothionamide**, **Rifampicin**, **Isoniazid**, and **Isoniazid High Dose**

Susceptible to: **Levofloxacin**, **Linezolid**, **Imipenem**, **Meropenem**, **Moxifloxacin**, **Moxifloxacin High Dose**, **Bedaquiline**, **Amikacin**, **Clofazimine**, **Terizidone**, **Streptomycin**, **Para Aminosalicylic Acid**, **Rifabutin**, **Delamanid**, and **Pyrazinamide**

Clinical Information

Hb (g/dL)	10.7	Platelets (*10 ⁹ /L)	600
Neutrophils (*10 ⁹ /L)	NA	GFR (mL/min)	60
ALT (IU/L)	13	QTc (ms)	391
Hearing loss	No	Painful peripheral neuropathy	No
Psychosis	No	Poorly controlled epilepsy	No
Visual problems	No	Pregnant	No

■ Clinical decision support system

- AI trained algorithm
- Takes DR-profile and clinical information into account
- Knowledge of drug properties

Treatment recommender:

Proposed Individual Treatment Regimen


Drug	Dosage (mg)
Levofloxacin	1000 mg
Linezolid	600 mg
Bedaquiline	200 mg three times per week
Pyrazinamide	1500 mg

Summary

- **NGS important role to accurately track transmission, outbreaks etc.**
- **Importance in WGS for drug resistance in TB:**
 - All known drug resistance conferring variants can be identified
 - Enables individualised treatment
 - Detection of heteroresistance
 - Identify novel mechanisms of resistance
- **Various applications of NGS**
 - Which clinical question are you seeking to address?

Tuberculosis practicals – Monday 2 Sep


- Full day!
- Please prepare – lab manual
- Before Thursday:
 - Create an account on Galaxy: <https://usegalaxy.org/>
 - Remember your login details...



COMPLEXITY

TORCH

Tuberculosis Omics Research



Stellenbosch
UNIVERSITY
IYUNIVESITHI
UNIVERSITEIT

TB Genomics
research group



torch-consortium.com



anzaan.dippenaar@uantwerpen.be

MOLECULAR ASSAYS FOR TUBERCULOSIS: increased levels of complexity

INFORMATION

Mtb/no *Mtb*
RIF

Mtb/no *Mtb*
RIF
INH

Mtb/no *Mtb*
15 drugs
Genotype

Mtb/no *Mtb*
all drugs
Genotype
Transmission

TARGETS

MTB-specific
1 drug-target

MTB-specific
3 drug-targets

MTB-specific
18 drug-targets

4.4Mb
All
drug-targets

SAMPLE
TYPE

Raw
specimen


Treated
specimen

Treated
Specimen

Cultured
isolate

MOLECULAR
ASSAY

Simple
NAATs



Moderate
complexity
NAATs

